

Comparison of PMA2100, DCS-1 and DCS2.0 in a SPF15 FDA Study using a 16S Single Port Simulator, increment 25%, from a base of 1 MED.

	PMA2100	DCS-1	DCS 2.0
Keystrokes: Set up	17-37	2-15	12
1 st test	3	9	2
2 nd test	3	7	2
3 rd test	3	5	2
4 th test	3	9	2
5 th test	3	6	2
Total (avg)	42	46	22
Daily (avg)	336	368	176
Ease of use:	Difficult	Easy	Easy